


BOOKWALTER™ STORAGE CASES

Now a better way to sterilize, store, and maintain your BOOKWALTER System and related components!



BOOKWALTER Table Fixation and Ring Case
Product Code #50-4507



BOOKWALTER Blades and Accessories Case
Product Code #50-4508



Codman

BOOKWALTER™ STORAGE CASES

<i>Features</i>	<i>Benefits</i>
Weight of system is evenly distributed	Less strain and fatigue to staff. Meets the Association for the Advancement of Medical Instrumentation (AAMI) guidelines.
Removable Cover	Easier access for operating room personnel
Ergonomically designed handle grips	Easier to carry
Laser etched product outlines & codes	Assists Central Sterile Department with assembly
High quality USA Stainless Steel	Corrosion resistant
Silicone rubber inserts	Protects products and facilitates sterilization process
Additional storage compartment	Provides storage for miscellaneous blades and components

ORDERING INFORMATION

50-4507 BOOKWALTER Table Fixation and Ring Case
Dimensions: 2-3/4" x 13-1/8" x 22-1/8"
Weight: 22-3/4 lbs.
Sterilization Method: Steam

These cases are designed to hold the following products:

50-4581 Table Post
50-4582 Horizontal Flex Bar
50-4552 Horizontal Bar
50-4583 Medium Segmented Ring
50-4557 Medium Oval Ring

50-4508 BOOKWALTER Blades and Accessories Case
Dimensions: 3-11/16" x 13-1/8" x 22-1/8"
Weight: 16-3/4 lbs.
Sterilization Method: Steam

These cases are designed to hold the following products:

50-4579 (2 ea.) Straight Ratchets
50-4580 (6 ea.) Tilt Ratchets
50-4554 Post Coupling
50-4566 Harrington Blade
50-4571 Kelly Blade 2x6
50-4570 Kelly Blade 2x5
50-4569 Kelly Blade 2x4
50-4568 (2 ea.) Kelly Blade 2x3
50-4576 Balfour Blade 3x2-5/8"
50-4577 Gelpi Blade
50-4567 Malleable Blade 3x6
50-4562 Malleable Blade 1-1/2x6
50-4561 (2 ea.) Malleable Blade 2x6

Codman
a Johnson & Johnson company

Raynham, MA 02767-0350 • (508) 880-8100
To Order Call: (800) 255-2500
For Product Information Call: (800) 225-0460

RET-10-000

06/06 QP/QP