


BOOKWALTER™ RETRACTOR ACCESSORIES

Important new additions to your CODMAN® BOOKWALTER Retractor Systems



The BOOKWALTER Retractor Systems, create an open field for the surgeon, enabling the operating room staff hands-free assistance

CODMAN BOOKWALTER RADIOLUCENT BLADE HOLDER

Benefits

- ❖ Holds the CODMAN BOOKWALTER Disposable Radiolucent Blades
- ❖ New Peak holder is radiolucent reducing artifacts from the holder during x-ray or fluoroscopy
- ❖ Compatible for use with all CODMAN BOOKWALTER components

Ordering Information

50-4677 CODMAN BOOKWALTER RADIOLUCENT BLADE HOLDER

Codman
a *Johnson & Johnson* company

See your Codman sales representative for information on the entire line of genuine BOOKWALTER Retractor products.

Raynham, MA 02767-0350 • 508-880-8100
To Order call: 800-255-2500
For Product Information Call: 800-225-0460
RET-25-000