

CODMAN[®] Precision Valve System


Proven Ball &
Spring Construction

Five Distinct Operating
Pressures


Multiple Configurations

Intergrated SIPHONGUARD[™]
Technology

Customised to Your Patients
Specific Needs



CODMAN Precision Valve long-term stability results. Numerous precision valves were submitted to a pulsatile flow test equivalent to four million activation cycles, demonstrating a high degree of long-term stability.*



Performance characteristics of competitive systems vs. the CODMAN Precision Valve. The CODMAN Precision Valve offers very refined operating bands without overlap.

*Data on file.

CODMAN® Precision Valve System

Valves

Inline with SIPHONGUARD

VL	825461
L	825462
ML	825463
MH	825464
H	825465



Inline

VL	825471
L	825472
ML	825473
MH	825474
H	825475



Right Angle with SIPHONGUARD

VL	825481
L	825482
ML	825483
MH	825484
H	825485



Right Angle

VL	825491
L	825492
ML	825493
MH	825494
H	825495



Cylindrical with Prechamber

VL	823011
L	823012
ML	823013
MH	823014
H	823015



Cylindrical

VL	823082
L	823083
ML	823084
MH	823085
H	823086



Micro Valve with RICKHAM® Reservoir

VL	823095
L	823096
ML	823097
MH	823098
H	823099



Micro Valve

VL	823026
L	823027
ML	823028
MH	823029
H	823030



Valve Systems

Includes: 14cm Ventricular Catheter, 120cm Unitized Distal Catheter

Inline with SIPHONGUARD

VL	823821
L	823822
ML	823823
MH	823824
H	823825



Inline

VL	823806
L	823807
ML	823808
MH	823809
H	823810



Right Angle with SIPHONGUARD

VL	823261
L	823262
ML	823263
MH	823264
H	823265



Right Angle

VL	823281
L	823282
ML	823283
MH	823284
H	823285



Includes: 14cm Ventricular Catheter, 85cm Unitized Distal Catheter

Cylindrical with Prechamber

VL	823006
L	823007
ML	823008
MH	823009
H	823010



Micro Valve with RICKHAM Reservoir

VL	823021
L	823022
ML	823023
MH	823024
H	823025



Micro Valve

VL	823035
L	823036
ML	823037
MH	823038
H	823039



Special Valves

Valves can be made to fit your special needs. Please contact your CODMAN Neuro Specialist or the CODMAN Specials Department at 1 (800) 843-0039 or log onto www.shuntworks.com.

Accessories: Catheters

823041	CODMAN Ventricular Catheter, 14cm
823044	CODMAN Atrial Catheter, 46cm
823045	CODMAN Peritoneal Catheter, 120cm
823072	BACTISEAL™ Catheter Kit
823073	BACTISEAL Ventricular Catheter, 14cm
823074	BACTISEAL Peritoneal Catheter, 120cm

Accessories: Miscellaneous

823048	CODMAN Stainless Steel Straight Connector
823049	CODMAN Stainless Steel Right Angle Connector
823052	CODMAN Plastic Right Angle Connector
823053	BIOPLATE™ Titanium Straight Connector
823055	CODMAN Plastic Valve Introducer
823090	SIPHONGUARD, CSF Flow Control Device



For more information, contact your Codman Sales Representative

325 Paramount Drive, Raynham, MA 02767 • (508) 880-8100
For product information, call (800) 225-0460 • www.codman.com

©2003 Codman & Shurtleff, Inc. All rights reserved.

VAL-14-001 8/03 DP/QP